

# 32" SmartLine Display Flexible Design and Size

## Unique Selling Points

- Adaptive frame design – suitable for various sizes of displays
- TrueFlat display solution - customized modifications possible
- Design available on long-term
- LED panel for brilliant colours and yet with low power consumption
- Full-HD 1920 × 1080
- High-intensity, high-contrast medical LCD
- Viewing angle – wide range of 178 degrees in both horizontal and vertical directions
- Multiple analog & digital I/O connectors




32"

### Short Reference



### 32" SmartLine Display



# 32" SmartLine Display Flexible Design and Size

## Technical Data

<b>Model / Type</b>	<b>32" SmartLine Display</b>
Screen size	32" / 80.01 cm
Active area (H × V)	698.4 × 392.85 mm
Resolution	1920 × 1080 / HDTV
Pixel pitch (H × W)	0.364 mm
Display colors	16.7 Mio. RGB
Brightness	350 cd/m <sup>2</sup> w/o touch
Contrast ratio	3000:1 typ.
Response time	6.5 ms typ.
Viewing angle: h / v	178° / 178° (CR=10)
I/O Ports	D-Sub 15pin (analog), DVI-D, DP on req: HDMI, Video & S-Video
Power consumption	max. 60 W
MTBF backlight	min. 30.000 h
Dimensions (W × H × D) w/o rotation arm adaption	758.4 × 457.4 × 80.1 mm
Weight	max. 15 kg
Operating temperature	0 °C to +40 °C
Storage temperature	-20 °C to +60 °C
Options	Protection shield, touch sensor, remote control, different colours, customer logo
Misc.	Speaker 2× 30 W, LED, slim, light weight, True Flat, protection shield, cabling accessible - but covered
Housing	Aluminium, colours, structure on request
Approval	CE, UL listed components
Mounting	VESA 200
Power Input	internal PS 24 V

## Different Views

